



Vibratory bowl not included.

StandPAK 5000 The perfect mate for Burr King Vibratory Bowls

- Accepts FilterPak 4001
- o-6 Hour Process Timer
  Timer Controlled AC
- outlets
- Isolator Stance
- Steel Construction Improves action
- Space Saving 16"x20"x21"
- Power Cord & Plug 110V, 60 Hz
- Mass finish many different material types, part shapes
- Perfect bench-top small mass finishing system
- Quiet operation
- Compatible with most media styles
- Compact footprint friendly to "cell" manufacturing

## VibraKING Bowl on the StandPAK 5000, add a FilterPak 4001 for a complete, timer controlled, filtered, compound finishing system.

Mount your

## VibraKING Vibratory Bowls & Accessories

SumpPAK 1000 comes with a 5 gallon bucket, pump hanger, pump and all tubes to connect filtration system to machine. Can be used with bowls or tubs.





The FilterPAK 4001 can be used with vibratory bowls or tubs. Heavy duty Polyethylene construction with single submersible pump and disposable filter bag. 10 gallon capacity. Also available FilterPAK 8000 used only with Vibra KING 85, stainless steel construction.

- Continuous fluid filtration
- Traps and holds sludge and debris in inexpensive, easily disposable filters
- Cleaner compounds results in cleaner parts
- Superior finishes and shorter cycle times with a filtered system
- Significantly reduces shop mess
- Filtration from 5 to 100 micron
- Filter change less than 2 minutes
- 115 VAC 140 GPH Pump
- 10 Gallon plus capacity
- Heavy duty Polyethylene
- Space saving 12"x14"x15"
- Dual Chamber

Always wear protective eye, face, arm, leg and body gear when using any grinder, polisher, buffer or other rotating machinery.BURR KING MANUFACTURING CO., INC.SINCE 19511220 TAMARA LANE,WARSAW, MO 65355(660) 438-8998(800) 621-2748FAX (660) 438-8991www.burrking.com

23